

350 overhead suspended shield & LED 150F lamp

350/45/LED150F-0031 (250 cm ceiling track, 40 cm x 50 cm lead glass shield & LED 150F operating lamp)

Track mounted, moveable upper body protection

This 40 cm x 50 cm shield uses 2.2 mm equivalent lead glass (@150 kV) for effective protection of the eyes and thyroid. The robust yoke assembly provides stability and protection to the window. The shield is fitted with a 0.5 mm Pb equivalent flexible curtain which is suspended beneath to further reduce scatter radiation. The LED 150F lamp is a small operating light with a light output of 110,000 lux.

350 x-ray shield, 40 cm x 50 cm

- > Protection for the eyes and thyroid.
- Landscape, 40 cm x 50 cm transparent shield with
 2.2 mm equivalent lead glass (@150 kV).
- Suspended, 0.5 mm Pb equivalent flexible curtain.
 20 cm in depth.
- > Pivots on the central horizontal axis.

LED 150F operating lamp

- > Light intensity at 1 metre: 110,000 Lux.
- > Colour rendering index R_a 97.
- > Colour rendering index R₉ 97.
- > Focusable light field size: 16 cm to 24 cm.
- Electronic light intensity control at lamp head: 50% to 100%.
- > Colour temperature: 4500 Kelvin.
- > Temperature increase at head area: 0.5°C.
- > Total power consumption: 24 W.
- > Light source: 26 LED's.
- > Life span of over 60,000 hours.
- > Diameter of lamp head: 40 cm.
- > Detachable sterilizable handle included.

Ceiling track system

- > Consisting of:
 - 250 cm aluminium track* and moveable carriage
 - 80 cm ceiling column*
 - 2 x two-part suspension arm
 - 40 x 50 cm shield
 - LED 150F lamp
- Optional extras:
 - Second track carriage for two additional suspended devices.
- Cable management system so lamp cables can be contained and managed within a separately supplied flexible energy chain.
- * See over other lengths available.




350/45/LED150F-0031

Kenex (Electro-Medical) Ltd Unit 17, RO24 Greenway, Harlow Business Park, Harlow, Essex CM19 5QB, England T: +44 (0)1279 417241 F: +44 (0)1279 443749 E: kenex@kenex.co.uk www.kenex.co.uk



Technical specification



Equipment weights & part numbers

350/45/LED150F-0031

Power supply in housing LED 150F lamp head: **3.5 kg** 40 cm x 50 cm lead glass shield: **12.5 kg** 2 x Suspension arm: **10.4 kg** (each) 80 cm dual mount ceiling column: **12.5 kg** 2.5 m ceiling track & carriage: **50 kg**

Optional equipment: Dual mount ceiling col

Dual mount ceiling column:	40 cm (#300/COLUMN/40): 10 kg
	60 cm (#300/COLUMN/60): 11.2 kg
	100 cm (#300/COLUMN/100): 14.7 kg
Ceiling track & carriage:	4.2 m (#3001/TC/420): 72 kg
Track cable management:	2.5 m (#3001/TC/250/EC): 1.3 kg
	4.2 m (#3001/TC/420/EC): 1.7 kg
Additional track carriage:	For use in 2.5 m or 4.2 m track (#3001/3.25): 4.6 kg

🚺 Kenex

Kenex (Electro-Medical) Ltd Unit 17, RO24 Greenway, Harlow Business Park, Harlow, Essex CM19 5QB, England T: +44 (0)1279 417241 F: +44 (0)1279 443749 E: kenex@kenex.co.uk www.kenex.co.uk